# Semen tank checklist

Don't let poor semen tank maintenance and management be the cause of poor semen fertility. Use this handy checklist to ensure proper handling for best semen performance.

Minimum level

should be

Use thermocoupler for optimum

temperature reading of

-220° F to -320° F

of liquid nitrogen

4 inches

## **Temperature Check**

#### WET TANK

Make a visual inspection of the tank. Check for excessive sweating and frost on outside of tank; excessive bubbling when removing neck cork. These could be signs of a potential issue.

Measure level of liquid nitrogen (N2). Level will be indicated by frostline, which develops when dipstick is removed.

#### **VAPOR TANK**

Make a visual inspection of the tank. Contact the shipping agent if any damage is apparent.

## **Inventory Management**

Create an inventory system to avoid searching for and potentially damaging or breaking straws.

### Tank Maintenance




NECK CORK: Should be a tight fit to avoid N2 leakage. A damaged or a loosefitting cork will result in N2 loss.

## Transfer

CANISTER: Stem or handle on the canister should not separate from the canister or bend.

the tank. Without it, there's a potential for tank failure.



CANISTER

Properly cool forceps or tweezers before handling of semen.

Tanks should be placed side-by-

**BLACK PLUG:** Needs to cover the valve at the back of

